



REITEL

# MINI-LAB

*Optimize your lab space with  
space-saving equipment combinations!*



- » Stable support through sturdy stainless steel casing
- » Efficient material removal with high-quality grinding discs
- » Easy-to-clean design with a removable front panel

☐ **ROTOGRIND CLASSIC** | 1.280 €  
Wet trimmer (without disc)

☐ **Set Corundum** | 1.500 €  
Trimmer + grinding disc

☐ **Set Sil.-Carb.** | 1.500 €  
Trimmer + grinding disc

☐ **Set PW** | 2.070 €  
Trimmer + grinding disc

☐ **Set DELTA** | 2.350 €  
Trimmer + grinding disc



- » Infinitely adjustable speed up to 4,500 RPM
- » Polishing unit can be wall-mounted with upward-directed air for full power
- » Easily replaceable and washable filter at 40°C

☐ **POLIRET MINI** | 2.700 €  
Polishing unit



- » Rapid pressure build-up with an integrated tank heater
- » Ergonomically shaped gun with a Luer Lock adapter
- » Versatile applications with various nozzle options



☐ **STEAMY MINI** | 1.450 €  
Steam cleaner

- » Optimal service for your patients with fast and effective cleaning performance
- » Infinitely adjustable speed control up to 3,000 RPM
- » The cup features an additional seal and is suitable for sterilization.

☐ **PROTOFRESH** | 1.060 €  
Denture cleaner



« See the experience report of a dental technician about his daily work with REITEL laboratory equipment.

- » Collects contaminants
- » Can be placed flexibly anywhere in the laboratory
- » Easy to clean thanks to a removable grid

☐ **STEAM BOWL** | 260 €  
Steam bowl



- » Cleans tools and instruments quickly and effectively – even in hard-to-reach areas
- » Temperature and time are adjustable individually
- » Comprehensive accessories available

☐ **SONIRET SWEEP** | 730 €  
Ultrasonic cleaning unit



« Please simply list the devices you are interested in and provide the company name along with the contact person.

- » No need for a compressed air connection, thanks to the integrated vacuum pump
- » Stirring with just one hand motion
- » Includes a 400 ml cup
- » Stand optional

☐ **VACURET EASY** | 1.450 €  
Vacuum mixing device




- » Independent laboratory workstation with ergonomically shaped real wood surface
- » Dust and gases are directly extracted
- » Includes vice and LED lighting
- » Optional magnifying glass

☐ **ERGORET COMFORT** | 950 €  
Workstation



☐ **WORK PAD** | 45 €  
Stainless steel



- » Three selectable light frequencies
- » Ample maneuverability with a large access surface
- » Timer can be set individually and to the second
- » Equipped with high-performance LEDs

☐ **POLYCURE** | 1.750 €  
Light curing device



- » Combines rotary and micro-fine blasting in a single device, saving space in the lab
- » Available with two or three abrasive chambers
- » Heated containers prevent nozzle clogging and ensure consistent material flow
- » Includes maintenance-free quick stop feature


☐ **SANDBLAST 2** | 2.400 €  
Recycling and microfine sandblaster



- » Brightness and color temperature are continuously adjustable
- » Perfect for accurate color reproduction at daylight level
- » Includes three memory slots for favorite light settings

☐ **LED MINI / COMFORT** | 890 € / 1.220 €  
Table installation set

☐ **SCREW CLAMP / WALL MOUNTING SET** | 100 € / 69 €  
Optional



- » Intensity can be adjusted individually
- » Two selectable vibration frequencies
- » Suction cup feet ensure stability and absorb vibrations
- » Includes a vibration/rattle plate

☐ **VIBROMAT** | 640 €  
Vibrator



Your purchase price\*

Your leasing rate\*

Distributor/Branch office\*

Address

Phone / Mail

Name / Signature