

■ Workstations

ERGORET CC

Workstation

Order-No.: 14803 (white) / 14803.1 (grey)



- ergonomic workstation, offering best conditions to minimise strain and tiredness
- potentially detrimental dust, gases and fumes can be safely extracted and are fed to the filter at the back of the unit
- additional extraction adapter for connection of an external extraction, such as the DUSTY COMFORT (page 8)
- massive wood base and protective glass cover
- including two drawers, trap for coarse waste, stainless steel tray, hinged magnifying glass, perfect LED lighting, air blast nozzle and adapter for handpiece
- function as a quality control / inspection unit
- optional Vario-lifting table (variable working height: 880 – 1.180 mm)

Technical data

Line voltage	230V / 50 – 60 Hz
Current consumption	120W
Height	670 mm
Width	910 mm
Depth	860 mm
Weight	60 kg
Magnifying glass	4 diopters
Compressed air	max. 8 bar
Housing	stainless steel / wood