

■ Cleaning

SUPERSTEAM

Steam cleaner

Order-No.: 14201



- dry steam, wet steam for cleaning and hot water spray jet for wax removal
- variable pressure, moisture and quantity of steam
- pressure 6 bar
- compressed air function to dry the object
- large tank
- pressure gauge
- easy filling and one-hand operation
- brushes for internal and external cleaning
- wall-mounted shelf is available as an optional extra

Technical data

Line voltage	230V / 50 – 60Hz
Current consumption	2.200 W
Height	420 mm
Width	300 mm
Depth	350 mm
Weight	14 kg
Temperature	160 °C
Pressure	6 bar
Capacity	5 l
Moisture + quantity	variable
Filling	manual filling
Housing	stainless steel