

# ■ Casting

## RETOCAST I

Induction casting centrifuge

**Order-No.: 18000**



- induction heating, variable performance
- optimally adapts to various alloys
- can be used for all common casting ring systems
- allows for quick and trouble-free casting
- easy handling, the intuitive control panel and the large sight glass offer a high degree of user comfort
- graphite or ceramic crucibles can be used
- coil moves out of the way when lid is opened
- electronic locked lid and timer
- emergency stop button
- three starting speeds

## Technical data

Line voltage	230V / 50 – 60Hz
Current consumption	3.200W
Height	530 mm open: 930 mm
Width	510 mm
Depth	740 mm
Weight	60 kg
Speed	500 rpm
Casting ring size	1 – 9 and model casting
Max. quantity of alloy	noble metal 110g non-precious metal 55g
Min. quantity of alloy	non-precious metal 15g
Housing	stainless steel