

■ Polishing

POLIRET TWIN

Polishing unit

Order-No.: 12703



- two brushless motors with adjustable speed control for continuous use
- POLIRET: fully variable speed
- POLIRET TWIN: two separate motors, each of them offering fully variable speed
- the polishing troughs made of plastic can easily be removed for cleaning
- easy change fleece filter, machine wash at 40°C
- emergency stop button with safety fuse

Technical data

Line voltage	230V/50–60 Hz
Motor power	2 x 370 W
Extraction power	400 W
Height	500 mm
Width	730 mm
Depth	560 mm
Weight	60 kg
Speed (variable)	800–4.500 rpm
Extraction capacity	800 m ³ /h
Housing	stainless steel