

# Mixing & Polymerisation



## POLYCURE

Polymerisation unit

**Order-No.: 18310000**



- three light frequencies for different materials
- durable high-performance LEDs
- wide access area and free movement thanks to hinged lid
- generously dimensioned rotary disc
- individually adjustable time with a digital timer
- integrated inner chamber ventilation prevents overheating of your workpieces

## Technical data

Line voltage	230 V / 50 – 60 Hz
Current consumption	100 W
Height	210 mm
	open: 440 mm
Width	210 mm
Depth	335 mm
Rotary disc-Ø	125 mm
Wavelength	385 – 390 nm
	405 – 410 nm
	465 – 470 nm
Weight	5,5 kg
Housing	stainless steel