

■ Mixing

VACURET CLASSIC

Vacuum mixing device

Order-No.: 10607100



- precise mixing of plaster, investment material and alginates with digital timer
- variable speed
- maintenance-free vacuum pump
- integrated, easy exchange filter system
- wall mounted, including 400 ml mixing bowl

Technical data

Line voltage	230V/50–60Hz
Current consumption	480 W
Height	250 mm
	incl. stand: 670 mm
Width	150 mm
	incl. stand: 215 mm
Depth	280 mm
	incl. stand: 355 mm
Weight	44kg
Speed	150 - 850rpm
Vacuum	max. 93 %
Sound level	≤ 70 dB(A)
Weight	8 kg / 15 kg inkl. stand
Housing	stainless steel