

Mixing

VACURET DIGIT

Vacuum mixing device with digital display

Order-No.: 10608100



- Vacuum mixer for plaster, investment material, silicone and certain resins, e.g. for 3D printing
- digital colour display with rotary knob
- maintenance-free vacuum pump
- brushless motor for continuous use
- variable speed 200-600 rpm, pre-mixing and main mixing adjustable
- 4 fixed programs, 28 programmable configurations, one manual setting
- incl. 400 ml plastic mixing bowl with stainless steel paddle and sturdy lid

Technical data

Line voltage	230V/50–60Hz
Current consumption	600 W
Height	250 mm
	incl. stand: 670 mm
Width	150 mm
	incl. stand: 215 mm
Depth	280 mm
	incl. stand: 355 mm
Weight	44kg
Speed	200 - 600rpm
Vacuum	max. 98 %
Sound level	≤ 70 dB(A)
Weight	8 kg / 15 kg inkl. stand
Housing	stainless steel