

# Mixing

## VACURET EASY INJECT

Vacuum mixing device

**Order-No.: 10606200**



- compressed air connection (without pump)
- easy mixing of plaster, investment material and alginates
- integrated easy exchange filter system
- wall mounted, including 400 ml mixing bowl

### Technical data

Line voltage	230V/50–60Hz
Current consumption	300 W
Height	250 mm
	<b>incl. stand: 670 mm</b>
Width	150 mm
	<b>incl. stand: 215 mm</b>
Depth	280 mm
	<b>incl. stand: 355 mm</b>
Weight	44kg
Speed	370 rpm
Vacuum	max. 90 %
Sound level	≤ 70 dB(A)
Weight	8 kg / 15 kg inkl. stand
Housing	stainless steel