

Mixing

RETOMIX MINI

Vacuum mixing device

Order-No.: 10607000



- for mixing plaster, investment material and alginates bubble-free and homogeneously
- electronic timer precise to the second
- integrated vacuum pump starts automatically
- variable vacuum intensity
- ease of operation
- easy-change filter system, preventing contamination and damage to the motor and the vacuum pump
- variable speed
- mixer is supplied for wall mounting
- incl. 400 ml mixing bowl
- mixing bowls are available in different sizes

Technical data

Current consumption	420 W
Line voltage	230 V / 50 – 60 Hz
Height	290 mm incl. stand: 710 mm
Width	155 mm incl. stand: 220 mm
Depth	285 mm incl. stand: 360 mm
Speed	variable 370 – 630 rpm
Vaakuum	max. 93 %
Weight	8 kg/15 kg mit Ständer
Housing	stainless steel