

Mixing & Polymerisation



VIBROMAT

Vibrator

Order-No.: 14002



- lateral vibration support
- designed to vibrate plaster, investment and other materials evenly and homogeneously
- with two variable vibration frequencies
- controllable intensity
- variable swinging magnet
- easy removable and washable rubber tray

Technical data

Line voltage	230 V / 50 – 60 Hz
Current consumption	150 W
Height	125 mm
Width	300 mm
Depth	170 mm
Weight	6 kg
Frequencies / min.	3.000 and 6.000
Intensity	continuously adjustable
Housing	stainless steel