

Cleaning

SONIRET RP HEAT

Ultrasonic cleaning unit

Order-No.: 10200006



- effective cleaning of tools and instruments with the new SWEEP Mode
- for optimum cleaning results in places which are difficult to access, e. g. bores, joints or grooves
- fast instrument cleaning
- high efficiency vibration technology due to automatic change of frequency
- integrated heating
- both time and temperature can be set on the control panel
- time-reduced disinfection procedure
- temperature range from 20 to 80 °C
- time setting adjustable from 1 – 99 minutes

Technical data

Line voltage	230 V / 50 – 60 Hz
Current consumption	650 W
Inhalt	2 x 5 l
<u>Device</u>	
Height	580 mm
Width	280 mm
Depth	297 mm
<u>Tank</u>	
Height	2 x 215 mm
Width	2 x 110 mm
Depth	2 x 215 mm
Weight	14,85 kg
Frequency	40 kHz
Housing	stainless steel