

Cleaning

SONIRET SWEEP

Ultrasonic cleaning unit

Order-No.: 10200100



- effective cleaning of tools and instruments with the new SWEEP Mode
- for optimum cleaning results in places which are difficult to access, e. g. bores, joints or grooves
- fast instrument cleaning
- high efficiency vibration technology due to automatic change of frequency
- integrated heating
- both time and temperature can be set on the control panel
- time-reduced disinfection procedure
- temperature range from 20 to 80 °C
- time setting adjustable from 1 – 99 minutes

Technical data

Line voltage	230 V / 50 – 60 Hz
Current consumption	200 W
Inhalt	max. 3 l
<u>Device</u>	
Height	210 mm
Width	280 mm
Depth	190 mm
<u>Tank</u>	
Height	100 mm
Width	240 mm
Depth	140 mm
Weight	4,5 kg
Frequency	39,5 – 40,5 kHz
Housing	stainless steel