

■ Tempering

THERMORET

Syringe / Composite warmer

Order-No.: 16002



- heats anesthetic and composite materials to allow warm delivery of the material as well as sodium hypochlorite solution and bite registration wax
- allows for easy composite flow
- material block enables preheating of carpules, syringes etc.
- easy-to-use display controller
- stainless steel water container

Technical data

Line voltage	230V / 50 – 60 Hz
Current consumption	180 W
Height	135 mm incl. box: 180 mm
Width	220 mm
Depth	230 mm
Weight	5,5 kg
Temperature	max. 70 °C
Capacity	approx. 0,8 l
Housing	stainless steel