

# ■ Trimming

## **ROTOGRIND MINI**

Wet trimmer

**Order-No.: 13400**



- precise finish with high surface quality
- adjustable table
- powerful brushless motor
- extremely stable and durable stainless steel housing
- front cover can easily be removed by two special quick releases for easy cleaning and fast changable discs
- safety fuse in the event of a power failure

### Technical data

Line voltage	230V / 50 – 60Hz
Current consumption	1.300W
Height	310 mm
Width	310 mm
Depth	400 mm
Weight	20 kg
Speed	3.000 rpm
Disc diameter	250 mm
Housing	stainless steel