

■ Trimming

ROTOGRIND NT KFO

Wet- / dry trimmer
with orthodontic angulation table

Order-No.: 43305



- perfectly suited for grinding precise angles onto maxillary and mandible models
- smooth axial movement
- high degree of freedom is provided for exact contouring
- if the table is not needed, it can easily be removed from the work area
- variable angle adjustment with locking positions
- extra large grinding disc opening (12 cm)
- model can also be precisely positioned using the laser guide (optional extra)
- comes with orthodontic front and grinding disc
- special diamond grinding disc for coarse and fine grinding - offering excellent surface quality

Technical data

Line voltage	230 V / 50 – 60 Hz
Current consumption	1.500 W
Height	350 mm
Width	350 mm
Depth	470 mm
Weight	22 kg
Speed	3.000 rpm
Disc diameter	280 mm
Housing	stainless steel