

■ Trimming

ROTOGRIND NT

Wet- / dry trimmer

Order-No.: 10301



- powerful wet trimmer, switchable to dry operation
- table can be fixed at any angle
- equipped with a silent brushless motor
- automatic motor stop and brake when opening the housing
- continuous water quantity control
- patented wet disc rinsing technology – even for models that are not fully hardened
- high stability – even for trimming super hard plasters
- hinged and removable front cover with quick release provides easy cleaning
- automatic motor stop and brake when opening the housing
- overload relay and safety fuse in the event of a power failure

Technical data

Line voltage	230 V / 50 – 60 Hz
Current consumption	1.500 W
Height	350 mm
Width	350 mm
Depth	470 mm
Weight	22 kg
Speed	3.000 rpm
Disc diameter	250 mm
Housing	stainless steel