


NEVER AGAIN BAD VIBRATIONS...

...with the new promotional set from REITEL!

Purchase SONIRET SWEEP by June 30, 2024, and receive a stainless steel insert basket for free as part of the set!



REITEL

SONIRET SWEEP WITH STAINLESS STEEL INSERT BASKET

Order-No.: 10200100
RRP: 730,- €



Details

Order-No.: 20400000
RRP: 42,- €



Ultrasonic cleaning is suitable for a wide range of materials.

Components that cannot withstand high temperatures can be effectively cleaned with ultrasonic technology. Even hard-to-reach areas are cleaned by the bubbles generated by the pressure and suction waves.

The SWEEP technique (frequency modulation) enables the creation of a particularly homogeneous and uniform ultrasonic field throughout the entire tank volume. The constantly changing frequencies prevent damage to delicate cleaning objects and ensure a highly effective and gentle cleaning process. The cleaning performance can be intensified as needed by activating the powerful integrated heater.

- » Temperature range from 20 to 80 °C
- » Switch-on time adjustable from 1 to 99 min.
- » Maximum tank capacity of 3 liters
- » New SWEEP technology with frequencies ranging from 39.5 to 40.5 kHz
- » Fast and effective cleaning even in hard-to-reach areas
- » Quick instrument turnover
- » Time-saving disinfection

Stainless steel basket:

- » Especially hygienic
- » Easy to hang
- » Robust and durable

Set price
730,- €
~~instead 772,- €~~

REITEL

REITEL Feinwerktechnik GmbH

Senfdamm 20 · 49152 Bad Essen · www.reitel.com · info@reitel.com