

# ■ Cleaning

## **BIOSTEAM CENTER**

Steam center

Order-No.: 10106



- clean and hygienic work environment
- the steam is absorbed at the backside of the steam center and condenses in the special filter
- absorbs potentially damaging air-borne by-products that are created when steam cleaning objects
- silent extraction
- bright LED lighting
- foot switch or manual switch for activation of the unit
- easy-change filter units
- drains on both sides usable
- installation independent of a basin

### Technical data

Line voltage	230 V / 50 – 60 Hz
Current consumption	200 W
Height	600 mm
Width	540 mm
Depth	510 mm
Weight	35 kg
Housing	Stainless steel