

# Mixing & Polymerisation



## **POLYMAT**

Polymerisation unit

**Order-No.: 18301**



- dental articulators or several flasks can be polymerised simultaneously
- pressure, time and curing temperature are adjustable
- heating process can also be started when the lid is open
- electronic temperature control
- timer scale up to 120 minutes
- heating outside the stainless steel tank
- status display
- fast pressure build-up
- water-proof according to IP 44
- various safety tools
- drain tap
- facilitated cleaning

## **Technical data**

Line voltage	230V/50–60 Hz
Current consumption	3.600W
Height	330 mm
Width	330 mm open: 290 mm
Depth	380 mm
Inner dimensions tank D/H	220/190 mm
Weight	22 kg
Capacity	7l
Pressure	0–6 bar
Temperature	max. 100°C
Housing	stainless steel