



■ DENTARY SURGERY EQUIPMENT



■ Cleaning

■ Tempering

ROTOCLEAN EASY

Brush cleaning for instruments

Order-No.: 10300



- suitable for manual pre-cleaning dental and surgical instruments prior to putting them into an autoclave
- direct contact with contaminated areas is completely avoided
- two built-in rotating brushes remove stubborn blood stains, secretions and calculus deposits
- brushes can be disinfected and sterilised up to 50 times
- protective glass against accidental splashing
- unit can be directly installed in the sink

SONIRET

Ultrasonic cleaning unit

Order-No.: 10200



- offers rapid and highly efficient precise cleaning method of different instruments and objects
- for optimum cleaning results in places which are difficult to access, e.g. bores, joints or grooves
- fast instrument cleaning
- high efficiency vibration technology
- integrated heating
- both time and temperature can be set on the control panel
- time-reduced disinfection procedure
- temperature range from 20 to 80 °C
- time setting adjustable from 1 – 99 minutes

THERMOMET

Syringe / Composite warmer

Order-No.: 16002



- heats anesthetic and composite materials to allow warm delivery of the material as well as sodium hypochlorite solution and bite registration wax
- allows for easy composite flow
- material block enables preheating of carpules, syringes etc.
- easy-to-use display controller
- stainless steel water container

Technical data

Line voltage	230V / 50 – 60Hz
Current consumption	30W
Height	160 mm open: 350 mm
Width	350 mm
Depth	270 mm open: 170 mm
Weight	6 kg
Capacity	2,7l
Housing	stainless steel

Line voltage	230V / 50 – 60Hz
Current consumption	200W
Inhalt	max. 3l
<u>Device</u>	
Height	210 mm
Width	280 mm
Depth	190 mm
<u>Tank</u>	
Height	100 mm
Width	240 mm
Depth	140 mm
Weight	4,5 kg
Frequency	40 kHz
Housing	stainless steel

Line voltage	230V / 50 – 60Hz
Current consumption	180W
Height	135 mm incl. box: 180 mm
Width	220 mm
Depth	230 mm
Weight	5,5 kg
Temperature	max. 70 °C
Capacity	approx. 0,8l
Housing	stainless steel

Recommendation

ROTOCLEAN EASY

Cleaning brush-set	20300
Hard cleaning brushes.....	20299
Cleaning brush on the front	20295

SONIRET

Insert tray, stainless steel.....	20400
Insert tray, plastic	20401
Beaker positioning cover, stainless steel	20403
Glass beaker	20402
Insert basket, plastic	20405
Insert basket, stainless steel	20404

OUR PERSONAL RECOMMENDATIONS

REITEL stainless steel care spray, 250 ml (Ingredients: animal oil, paraffin oil, inhibitors).....	90405
-------------------------------------------------------------------------------------------------------	-------